

APPLICATION NOTE

Monitoring Triazines and their Precursors Using Fiber Optic-based, UV-VIS Spectroscopy

Our product line includes a 508 UV-VIS™ spectrometer. This application note discusses the use of our hardware and software tools for the measurement of various triazines using fiber optic-based, UV-VIS spectroscopy. UV-VIS spectroscopy can be applied in real time directly in process monitoring or as a laboratory procedure. In either case UV-VIS spectroscopy is a time and money saving alternative to traditional chemical methods.




Figure 1: Absorbance Spectra of Atrazine Baseline Corrected at 330 and 530nm


Background

Atrazine, a triazine, is a restricted use pesticide in the United States. It is one of the most commonly used herbicides in the US particularly in the mid-western “corn belt”. It is considered slightly toxic, believed to be an endocrine disruptor and is a suspect carcinogen. Atrazine is a common contaminant in ground and surface water where it is slow to degrade. Agricultural runoff is the major source of surface water contamination. The US EPA Maximum Containment Level Goal (MCLG) for drinking water is 3ppb.

The UV-VIS region of the electromagnetic spectrum displays electronic transitions and is particularly useful for viewing conjugated and aromatic molecules. By measuring the UV spectra of a series of samples of known triazine concentrations, quantitative models can be developed which will allow the measurement of future samples based only on their UV spectrum. Our analyzer systems use fiber optics to allow the sample probe to be located in remote locations away from the spectrometer itself, potentially reducing the level of operator intervention and providing real time sample assessment.

Experimental

The triazines were initially diluted in isopropanol and then serially diluted with deionized water. The UV spectra were measured using a 508 UV-VIS process spectrometer in a 100mm path length sample cell. Isopropanol is UV transparent. The unit was referenced using water, allowing the spectral features attributed to the atrazines to be more easily seen. Figure 1 shows the absorbance spectra collected for atrazine itself baseline corrected at 330 and 350nm.



Model 508 UV-VIS Analyzer

Analysis

First derivative spectra preprocessing was used to eliminate baseline shifts. PLS models were then created for each chemical. Results shown in table below.

Chemical	Correlation	RMSEP wt/wt ppb
Atrazine	.995	9.11
Propazine	.999	2.57
Cyanuric Chloride	.996	11.35




Figure 2: Absorbance Spectra of Propazine baseline corrected at 330 and 530 nm

Discussion

The measurement of ppb level atrazine, propazine, or cyanuric chloride concentrations in water using UV spectroscopy is both fast and reliable utilizing the 508 UV-VIS analyzer. This method minimizes the need for laboratory sample collection. Results are available in real-time (seconds) for triazine concentrations in aqueous streams. For more detailed information regarding system specifications please contact a Process Insights technical sales specialist.




Figure 3: Absorbance Spectra of Cyanuric Chloride, a Triazine Precursor, Baseline Corrected at 330 and 530 nm

GAIN REAL-TIME INSIGHT INTO YOUR PROCESS

Process Insights manufactures and delivers premium sensors, monitors, detectors, analyzers, instrumentation, and software that are mission-critical to keep your operations, personnel, and the environment safe – every day across the globe.

Get the most reliable, precision analytical technologies available on the market today. We will work to match your needs and budget, and provide the optimal, and most stable process analysis solution for your application.

CENTERS OF EXCELLENCE | PROVIDING PROVEN SOLUTIONS

Process Insights is committed to solving our customers' most complex analytical, process, and measurement challenges everyday.

Process Insights – The Americas

4140 World Houston Parkway Suite 180, Houston, TX 77032, USA +1 713 947 9591

Process Insights – EMEA

ATRICOM, Lyoner Strasse 15, 60528 Frankfurt, Germany +49 69 20436910


Process Insights – APAC

Wujiang Economic and Technology, Development Zone, No. 258 Yi He Road, 215200 Suzhou, Jiangsu Province, China +86 400 086 0106

For a complete range of products, applications, systems, and service options, please contact us at: info@process-insights.com

For a complete list of sales & manufacturing sites, please visit: <https://www.process-insights.com/about-us/locations/>

COSA Xentaur, Tiger Optics, Extrel, Alpha Omega Instruments, ATOM Instrument, MBW Calibration, MGA, Guided Wave, ANALECT and LAR TOC Leader are trademarks of Process Insights, Inc.



REVOLUTIONIZING MEASUREMENT

EVERYWHERE